

<u>Darby sku</u>	<u>Description</u>	<u>SHM</u>	<u>Color</u>	<u>Work Time</u>	<u>Intra-Oral Set Time</u>	<u>Total Time</u>	<u>Shore "A" Durometer</u>
5252801	Alginate Substitute Fast Set, 50 ml	Yes	Pink	1:00	1:00	2:00	45
5252802	Alginate Substitute Fast Set, 380 ml	Yes	Pink	1:30	3:00	4:30	45
5252799	Alginate Substitute Regular Set, 50 ml	Yes	Pink	2:30	1:00	3:30	45
5252800	Alginate Substitute Regular Set, 380 ml	Yes	Pink	2:00	3:30	5:30	45
5252754-56	Bite Registration Super Fast Set	No	Light Blue	0:30	0:45	1:15	85
5252750	Bite Registration Fast Set, Rigid	No	Light Blue	0:40	1:00	1:40	93
5252751-53	Bite Registration Fast Set	No	Light Blue	0:40	1:00	1:40	85
5252748	Bite Registration Regular Set, Clear	No	Clear	1:00	2:00	3:00	70
5252749	Bite Registration Regular Set, Unflavored	No	Light Blue	2:00	1:00	3:00	85
5252761-64	Light Body, Fast Set	Yes	Yellow	1:00	1:00	2:00	50
5252769-72	Medium Body, Fast Set	Yes	Tan	1:00	1:00	2:00	55
5252778-79 / 81-82	Heavy Body, Fast Set	Yes	Rose	1:00	1:00	2:00	75
5252780	Heavy Body Fast Set, 380 ml Refill	Yes	Rose	1:30	3:00	4:30	75
5252788-89 / 91-92	Monophase, Fast Set	Yes	Purple	1:00	1:00	2:00	60
5252790	Monophase, Fast Set, 380 ml Refill	Yes	Purple	1:30	3:00	4:30	60
5252757-60	Light Body, Regular Set	Yes	Yellow	2:30	1:00	3:30	50
5252765-68	Medium Body, Regular Set	Yes	Tan	2:30	1:00	3:30	55
5252773-74 / 76-77	Heavy Body, Regular Set	Yes	Rose	2:30	1:00	3:30	75
5252775	Heavy Body, Regular Set, 380 ml Refill	Yes	Rose	2:00	3:30	5:30	75
5252783-84 / 85-86	Monophase, Regular Set	Yes	Purple	2:30	1:00	3:30	60
5252785	Monophase, Regular Set, 380 ml Refill	Yes	Purple	2:00	3:30	5:30	60
5252796-98	Putty, Fast Set	No	Dark Blue	1:00	1:00	2:00	72
5252793-95	Putty, Regular Set	No	Dark Blue	2:30	1:00	3:30	72