



QD NON-OPTIC COUPLERS



Darby 4-Hole
9546084



Beyes 4-Hole
9558810



Vector 4-Hole
***Special Order**

Compatible High Speed Handpieces:

9546083 - Darby Super Trac Pro Mini Head Push Button

9550039 - Beyes AirLight M800 Plus

Compatible High Speed Surgical Handpiece:

9490054 - Vector ECO45 Surgical High Speed (Includes QD Coupler)

BEYES PD NON-OPTIC COUPLER



Beyes 4-Hole
9558833

Compatible High Speed Handpieces:

9550033 - Beyes AirLight M800 Plus

9550095 - Beyes AirLight M800X

5251283 - Beyes DuraMax M800-NS

Compatible High Speed Surgical Handpiece:

9550054 - Beyes AirLight M800 Surgical 45

To order, call **800.645.2310** or visit **darby.com**