



KAVO/KAVO STYLE MULTIFLEX FIBER OPTIC COUPLERS



KaVo 5-Hole
8700614

KaVo 6-Pin
LED
8701030

Darby
5-Hole
9546468

Darby 6-Pin
LED
9546469

NSK 6-Pin
LED
9543138

Dentsply
5-Hole
8642472

Dentsply
6-Pin LED
8642474

Beyes 6-Pin
LED
9558836

Vector 6-Pin
9490060

Compatible High Speed Handpieces:

- 9546466** - Darby Super Trac Prestige Fiber Optic High Speed
- 8700711** - KaVo MASTERtorque M8900 L
- 8700775** - KaVo MASTERtorque M8700 L
- 8701041** - KaVo EXPERTtorque E679L
- 8701040** - KaVo EXPERTtorque E675L
- 8701046** - KaVo SMARTtorque S619L
- 8701044** - KaVo SMARTtorque S615L
- 5255664** - NSK Ti-Max Z890KL
- 5255667** - NSK Ti-Max Z990KL
- 5255662** - NSK Ti-Max Z Micro KL
- 9543144** - NSK Ti-Max X600KL
- 9543201** - NSK Ti-Max X500KL

- 9543224** - NSK S-Max M600KL
- 9543223** - NSK S-Max M500KL
- 8640357** - Dentsply Midwest Phoenix ZR
- 8640354** - Dentsply Midwest Phoenix
- 8640468** - Dentsply Midwest Stylus Plus
- 8642278** - Dentsply Midwest Tradition Pro
- 8942210** - Star 430SWL Torque Flex Lube Free
- 8942211** - Star 430SWL Torque Flex Lube Needed
- 9551917** - Morita TwinPower Turbine Standard Head
- 9490020** - Vector F6-SLK Torque Head
- 9490010** - Vector F4-SLK Ultra Access Head
- 5250139** - Beyes Maxso X200

Compatible High Speed Surgical Handpieces:

- 8700134** - KaVo SURGtorque S459L
- 9543148** - NSK Ti-Max X450KL
- 9551665** - Vista Impact-Air 45 Fiber Optic

To order, call **800.645.2310** or visit **darby.com**



NSK FIBER OPTIC COUPLERS



NSK 5-Hole
9543135



NSK 6-Pin LED
9543137

Compatible High Speed Handpieces:

5255663 - NSK Ti-Max Z890L

5255666 - NSK Ti-Max Z990L

5255661 - NSK Ti-Max Z Micro

9543143 - NSK Ti-Max X600L

9543199 - NSK Ti-Max X500L

9543222 - NSK S-Max M600L

9543221 - NSK S-Max M500L

8640353 - Dentsply Sirona Midwest Phoenix

8640467 - Dentsply Sirona Midwest Stylus Plus

9559740 - Morita TwinPower Turbine Standard Head

5250140 - Beyes Maxso X200

Compatible High Speed Surgical Handpiece:

9543147 - NSK Ti-Max X450L

NSK NON-OPTIC COUPLER



NSK 4-Hole
9543136

Compatible High Speed Handpieces:

9543226 - NSK S-Max M600

9543225 - NSK S-Max M500

9550096 - Beyes Airlight M800X

9550034 - Beyes Airlight M800 Plus

Compatible High Speed Surgical Handpieces:

9543145 - NSK Ti-Max X450

9550053 - Beyes AirLight M800 Surgical 45

To order, call **800.645.2310** or visit **darby.com**



DENTSPLY MIDWEST FIBER OPTIC COUPLERS



Dentsply Midwest 5-Hole
8642510



Dentsply Midwest 6-Pin
8642512

Compatible High Speed Handpieces:

8640356 - Dentsply Sirona Midwest Phoenix ZR

8640351 - Dentsply Sirona Midwest Phoenix

8640465 - Dentsply Sirona Midwest Stylus Plus

To order, call **800.645.2310** or visit **darby.com**



DENTSPLY SIRONA FIBER OPTIC COUPLER



Dentsply Sirona 5-Hole
***Special Order**

Compatible High Speed Handpieces:

8640352 - Dentsply Sirona Midwest Phoenix

8640466 - Dentsply Sirona Midwest Stylus Plus

DENTSPLY SIRONA NON-OPTIC COUPLER



Dentsply Sirona 4-Hole
***Special Order**

Compatible High Speed Handpiece:

8642275 - Dentsply Sirona Midwest Tradition TC

To order, call **800.645.2310** or visit **darby.com**