



# KAVO/KAVO STYLE MULTIFLEX FIBER OPTIC COUPLERS



KaVo 5-Hole  
8700614



KaVo 6-Pin  
LED  
8701030



Darby  
5-Hole  
9546468



Darby 6-Pin  
LED  
9546469



NSK 6-Pin  
LED  
9543138



Dentsply  
5-Hole  
8642472



Dentsply  
6-Pin LED  
8642474



Beyes 6-Pin  
LED  
9558836



Vector 6-Pin  
9490060

## Compatible High Speed Handpieces:

- 9546466** - Darby Super Trac Prestige Fiber Optic High Speed
- 8700711** - KaVo MASTERtorque M8900 L
- 8700775** - KaVo MASTERtorque M8700 L
- 8701041** - KaVo EXPERTtorque E679L
- 8701040** - KaVo EXPERTtorque E675L
- 8701046** - KaVo SMARTtorque S619L
- 8701044** - KaVo SMARTtorque S615L
- 5255664** - NSK Ti-Max Z890KL
- 5255667** - NSK Ti-Max Z990KL
- 5255662** - NSK Ti-Max Z Micro KL
- 9543144** - NSK Ti-Max X600KL
- 9543201** - NSK Ti-Max X500KL

- 9543224** - NSK S-Max M600KL
- 9543223** - NSK S-Max M500KL
- 8640357** - Dentsply Midwest Phoenix ZR
- 8640354** - Dentsply Midwest Phoenix
- 8640468** - Dentsply Midwest Stylus Plus
- 8642278** - Dentsply Midwest Tradition Pro
- 8942210** - Star 430SWL Torque Flex Lube Free
- 8942211** - Star 430SWL Torque Flex Lube Needed
- 9551917** - Morita TwinPower Turbine Standard Head
- 9490020** - Vector F6-SLK Torque Head
- 9490010** - Vector F4-SLK Ultra Access Head
- 5250139** - Beyes Maxso X200

## Compatible High Speed Surgical Handpieces:

- 8700134** - KaVo SURGtorque S459L
- 9543148** - NSK Ti-Max X450KL
- 9551665** - Vista Impact-Air 45 Fiber Optic

To order, call **800.645.2310** or visit **darby.com**



## NSK FIBER OPTIC COUPLERS



NSK 5-Hole  
9543135



NSK 6-Pin LED  
9543137

### Compatible High Speed Handpieces:

**5255663** - NSK Ti-Max Z890L

**5255666** - NSK Ti-Max Z990L

**5255661** - NSK Ti-Max Z Micro

**9543143** - NSK Ti-Max X600L

**9543199** - NSK Ti-Max X500L

**9543222** - NSK S-Max M600L

**9543221** - NSK S-Max M500L

**8640353** - Dentsply Sirona Midwest Phoenix

**8640467** - Dentsply Sirona Midwest Stylus Plus

**9559740** - Morita TwinPower Turbine Standard Head

**5250140** - Beyes Maxso X200

### Compatible High Speed Surgical Handpiece:

**9543147** - NSK Ti-Max X450L

## NSK NON-OPTIC COUPLER



NSK 4-Hole  
9543136

### Compatible High Speed Handpieces:

**9543226** - NSK S-Max M600

**9543225** - NSK S-Max M500

**9550096** - Beyes Airlight M800X

**9550034** - Beyes Airlight M800 Plus

### Compatible High Speed Surgical Handpieces:

**9543145** - NSK Ti-Max X450

**9550053** - Beyes AirLight M800 Surgical 45

To order, call **800.645.2310** or visit **darby.com**



## W&H RQ-24 FIBER OPTIC COUPLER



W&H 6-Pin  
**8701519**

### Compatible High Speed Handpieces:

- 8701501** - W&H Synea Vision TK-97L
- 8701503** - W&H Synea Vision TK-98L
- 8701500** - W&H Synea Fusion TG-97L
- 8701502** - W&H Synea Fusion TG-98L
- 5255665** - NSK Ti-Max Z890WL
- 5255668** - NSK Ti-Max Z990WL

## W&H RQ-04 NON-OPTIC COUPLER



W&H 4-Hole  
**8701518**

### Compatible High Speed Handpieces:

- 5255761** - W&H Synea Fusion TG-97
- 5255762** - W&H Synea Fusion TG-98
- 9550040** - Beyes AirLight M800 Plus

To order, call **800.645.2310** or visit **darby.com**