




# DARBY PARALLEL FIBER POSTS



Darby Parallel Fiber Posts are glass fiber posts that support restorations and prosthetic crowns. They feature a parallel shape with grooves for great retention and metal filament that allows for radiographic visualization.

**9520867** Starter Kit, 20 Parallel Posts, 5 each size and 1 drill of each size

Intro Kits: 10 Parallel Posts and 1 Drill

**9520868** Size 1 Kit, 0.9mm, White

**9520869** Size 2 Kit, 1.1mm, Yellow

**9520870** Size 3 Kit, 1.3mm, Red

**9520871** Size 4 Kit, 1.5mm, Blue

Refill Kits: 10 Parallel Posts

**9520872** Size 1 Refill, 0.9mm, White

**9520873** Size 2 Refill, 1.1mm, Yellow

**9520874** Size 3 Refill, 1.3mm, Red

**9520875** Size 4 Refill, 1.5mm, Blue

## Features and Benefits:

- Retentive shape
- Longitudinal fibers for easy removal
- High concentration of fibers and great mechanical properties
- Modulus of elasticity similar to dentin lowers risk of root fracture
- Conical tip for less tissue removal in apical third of root canal
- Radiopaque

To order, call **800.645.2310** or visit **darby.com**